

A large, thick, orange brushstroke graphic that curves from the top right towards the bottom left, framing the title text.

# Atertia® SafeControl® Distributor Pack

**Atertia® SafeControl®**

# Contents

---

The Ateria® SafeControl® vision	3
Market segmentation 1 – Standard pen needle users	4
Market segmentation 2 – Passive pen needle users	5
Product overview	6
Device specifications	8
Packaging	9
Packing and shipping configurations	10
Instructions for use	11
Compatibility	12
Order codes	13
References	14

## **The Ateria® SafeControl® vision**

### **A preference over conventional and safety pen needles**

We aim to establish Ateria® SafeControl® as the next generation of safety pen needles that is preferred by Healthcare Professionals and buyers because it provides a balance of safety and control during the injection process.

### **Unique selling proposition (USP)**

Ateria® SafeControl® is the next generation of safety pen needles. It is designed to protect Healthcare Professionals from needlestick injuries and provide confidence in delivering a full medication dose by balancing both safety and control during the injection process.

## Market segmentation 1

Standard pen needle users – Objective: to convert from standard pen needles to safety pen needles.

### Main challenge of standard pen needles

Healthcare Professionals (HCPs) are at risk of needlestick injuries.<sup>1</sup>



Needlestick injuries can have potentially serious physical and psychological consequences

### Solution: Features and benefits of using Ateria® SafeControl®

1. Ateria® SafeControl® is a safety pen needle so protects from needlestick injuries.



**Safety guard** for covering the needle after injection to protect from needlestick injuries

2. Ateria® SafeControl® uses the same injection technique as conventional pen needles.



**Visible needle** for control during the injection process and using the same technique as with conventional pen needles

## Market segmentation 2

Passive safety pen needle users – objective; switching from passive to manual activation.

### Challenges of passive safety pen needles

Limited control during the injection process as users are uncertain when the passive safety mechanism activates.<sup>2</sup>



Limited control during the injection process<sup>2</sup>

Premature activation of the safety mechanism resulting in uncertainty that the full medication dose has been delivered to the patient.



Premature activation of the safety mechanism<sup>2</sup>

Can't see the needle entering the patient's skin leading to uncertainty of the dose delivery.



Lack of needle visibility<sup>2</sup>

### Solution: Features and benefits of using Ateria® SafeControl®

1. Innovative manual push tab means Healthcare Professionals are both protected from needlestick injuries and are in control of the injection process resulting in confidence the full dose has been delivered.



**Innovative manual push tab** for easy, one-handed activation of the safety mechanism

2. Visible needle allows Healthcare Professionals see the needle going into the skin, resulting in confidence the full dose has been delivered.



**Visible needle** for control during the injection process and using the same technique as with conventional pen needles

## Product overview

The next generation of safety pen needles, designed to protect Healthcare Professionals from needlestick injuries and provide confidence in delivering a full medication dose by balancing both safety and control during the injection process.



# Ateria® SafeControl® Support

- We have a full suite of pre-launch and launch marketing communication collateral available for your market.
- These materials have been developed in conjunction with global market research<sup>2</sup> supporting Ateria® SafeControl® product claims.
- Please speak to your local Owen Mumford representative for information.



Ateria® SafeControl® brochure



Ateria® SafeControl® whitepaper

## Device specifications



Product data (legislation and reference standards)	
<b>CE classification</b>	Ila
<b>Sterilisation</b>	Gamma irradiation
<b>Validity</b>	5 years
<b>Biocompatibility</b>	ISO 10993 (2009)
<b>Legislation and reference standards</b>	ISO 11608-2, ISO 11137, ISO 11607

Pen needle components and relative raw materials	
<b>Hub (attaches to pen)</b>	Polypropylene
<b>Safety guard</b>	High impact polystyrene
<b>Dual safety indicator</b>	Polypropylene
<b>Green cap</b>	Polypropylene
<b>Cannula (needle)</b>	Stainless steel AISI 304

Dimensional characteristics				
Needle length (mm)	Gauge	Cannula external diameter (mm)	Cannula internal diameter (mm)	Wall thickness (ISO 9626-2)
5	30G	0.30	0.19	Extra thin wall
8	30G	0.30	0.19	Extra thin wall



## Packaging

There are seven language sets of packaging per variant. Please refer to page 13 or speak to your local Owen Mumford representative for information.




5mm packaging




8mm packaging

## Packing and shipping configuration




### Cartons

<b>Total number of devices per carton</b>	100	
<b>Device weight</b>	4.3g per unit	
<b>Carton dimensions</b>	Width	139.5mm
	Height	149mm
	Depth	111.75mm / 111.25mm
<b>Total weight of one individual carton</b>	430g	



### Shipping cases

<b>Configuration of a case</b>	20 cartons in a shipping case	
<b>Total number of devices in a case</b>	100 per carton x 20 cartons = 2,000 per case	
<b>Shipping case dimensions</b>	Width	579mm
	Height	284mm
	Depth	299mm
<b>Total weight of a filled case</b>	2,000 devices x 4.3g/device = 8.6kg	




### Shipping pallets

<b>Configuration of a pallet</b>	24 cases on a pallet, 4 layers of 6 cases	
<b>Total number of devices on a pallet</b>	2,000 per case x 24 cases = 48,000 per pallet	

## Instructions for use

1. Remove protective seal. Screw needle onto injection pen and remove green cap. Inject as per pen instructions.
2. After injection, push the white safety guard away from you until it clicks and covers the needle. You will now see an orange indicator.
3. Unscrew needle from the injection pen and dispose of in sharps container.



How-to-use video can be seen at: [safecontrol.com](http://safecontrol.com)

## Compatibility

[owenmumford.com/compatibility](http://owenmumford.com/compatibility)



Including:



Autopen® Classic



Autopen® 24



Including:



Byetta®



Including:



Forsteo® / Forteo®



KwikPen™



HumaPen® Luxura™ HD



HumaPen® Savvio™



Including:



Fiasp® FlexTouch®



NovoPen® FlexPen®



NovoPen® FlexTouch®



NovoPen Echo®



NovoPen® 5



Victoza®



Including:



AllStar™



Apidra® SoloStar®



ClikSTAR®



Lantus® SoloStar®



Toujeo® SoloStar®

Compatibility data correct as of March 2019.

## Order codes and language sets

Product description	Code	Devices per carton	Cartons per case	Devices per case	Case weight (kg)	Case size (cm)	Pallet config	Devices per pallet
<b>5mm x 0.30mm (30G)</b>	See table below	100	20	2000	8.6	<b>Width:</b> 579mm	24 (4 layers of 6 cases)	48,000
<b>8mm x 0.30mm (30G)</b>						<b>Height:</b> 284mm <b>Depth:</b> 299mm		

Needle length mm	Needle gauge	Description	Language set	Product code
<b>5</b>	30	Aterial® SafeControl® 30G 5mm CE	Bulk	001870
			1A	002039
			1B	002040
			2A	002041
			2B	002042
			3A	002043
			3B	002044
			4 (Malaysian)	002045
<b>8</b>	30	Aterial® SafeControl® 30G 8mm CE	Bulk	001869
			1A	002046
			1B	002047
			2A	002048
			2B	002049
			3A	002050
			3B	002051
			4 (Malaysian)	002052

<b>Set 1A</b>	English (EN)	French (FR)	N Netherlands (NL) (Dutch)	German (DE)	Italy (IT) (Italian)	Portugal (PT) (Portuguese)
<b>Set 1B</b>	English (EN)	Denmark (DA) (Danish)	Norway (NO) (Norwegian)	Swedish (SV)	Finland (FI) (Finnish)	Spain (ES)
<b>Set 2A</b>	English (EN)	French (FR)	Spain (ES)	Portugal (PT) (Portuguese)	Russia (RU) (Russian)	Arabic (AR)
<b>Set 2B</b>	English (EN)	Hebrew (HE)	Greek (EL)	Turkey (TR) (Turkish)	Ukraine (UK) (Ukrainian)	
<b>Set 3A</b>	Bulgaria (BG) (Bulgarian)	Croatia (HR) (Croatian)	Czech (CS) (Czech)	Estonia (ET) (Estonian)	Hungary (HU) (Hungarian)	Latvia (LV) (Latvian)
<b>Set 3B</b>	Lithuania (LT) (Lithuanian)	Romania (RO) (Romanian)	Slovenian (SL)	Slovakia (SK) (Slovakian)	Serbian (SR)	Poland (PL) (Polish)
<b>Malaysian</b>	English (EN)	Malaysia (MY)	Indonesia (ID) (Indonesian)	China (CN) (Simplified Chinese)	Thailand (TH) (Thai)	Vietnam (VN)

## Ordering information

### Contact details:

#### Asia Pacific

Shirley Loh – Regional Managing Director

T : +60 (0)7 553 9617

E : shirley.loh@owenmumford.com

#### DACH – Germany, Austria and Switzerland

Thomas Roßbach – Regional General Manager

T : +49 6026 97750

E : thomas.rossbach@owenmumford.com

#### North America and Caribbean

Travis Shaw – Executive VP and General Manager

T : +1 770 977 2226

E : t.shaw@owenmumfordinc.com

#### Northern Europe and Rest of World

Joanne Thompson – Regional General Manager

T : +44 1608 645 555

E : joanne.thompson@owenmumford.com

#### Southern Europe & Northern Africa

Eric Del Cotto – Regional General Manager

T : +33 232 518 870

E : eric.delcotto@owenmumford-fr.com

### Ordering lead times:

Please contact your local Owen Mumford representative for information.

### Pricing:

Please contact your local Owen Mumford representative for information.

#### References

1. Motaarefi H, Mahmoudi H, Mohammadi E, Hasanpour-Dehkordi A (2016) Factors Associated with Needlestick Injuries in Health Care Occupations: A Systematic Review. *J Clin Diagn Res* 10:1E01–1E04. 2. Project Saturn A (2017) Online study commissioned with an independent market research agency. Data on file.