(en)

Indications: For use on adhesives, dressings and appliances on skin.
DESCRIPTION OF DEVICE & CONTENTS: ESENTA™ Medical adhesive remover spray sting free (50ml)

INTENDED USE: Removal of adhesives, dressings and appliances from skin.

INSTRUCTION FOR USE:

- Spray around the edge of the adhesive and wait a few seconds.
- Gently peel adhesive away from the skin using additional spray to ease removal if necessary.

INFORMATION: Dries in a few seconds. Sprays at any angle. If during the use of this device or as a result of its use a serious incident has occurred, please report it to the manufacturer and to your national authority.

CAUTION:

- Avoid direct use on wound beds.
 - Avoid use on delicate or sensitive areas of skin (eyes, mouth, ears, nose, vagina).
 - If in doubt, patch test a small area prior to use.
 - Should skin irritation occur, discontinue use and consult a healthcare professional.
 - Avoid inhalation.
- For external use only.

WARNINGS: Use only as directed. Pressurised container: may burst if heated. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Do not pierce or burn even after use. Protect from sunlight. Do not expose to temperatures exceeding 50°C/122°F. Keep out of reach of children. EXTREMELY FLAMMABLE

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ARTWORK SPECIFICATION	
Artcode	SI2389
Artwork Description	INS ESENTA ADHES REMOVER SPRAY(50ML)INT
ICC	423289
SAP	1729206
KOM/REF	N/A
IS Code	N/A
Component Type	Insert
LR#	LR-004247
Existing Art Code being Revised	N/A
REQUIREMENT	SPECIFICATION
Bucket	INT
Local Code(s)	N/A
Size Code / Dieline	N/A
Dimensions	Flat: 363 mm x 200mm, Folded: 33mm x 200mm
Color(s) / Colour(s)	Black
Specification Reference(s)	TP/TF21-054
Bar Code: UPC (A) or GTIN*	N/A
Bar Code: PZN – Code 39*	N/A
Bar Code: JAN-13*	N/A
* NOTE: Bar Code: See artwork for positioning, magnification and truncation. Scanability verified to ANSI grade “C” or better.	
HISTORY	
28-Jan-2021 / SI2389: New insert for Esenta Adhesive Remover spray.	

Workflow Name	Workflow Start Date	User	Intention	Date	Justif.
Workflow - LBL Approval with Proof Review Apr 7, 2021	07-Apr-2021 05:51:40 (GMT)	Iveta Bajůžiková	Quality Approval	11-May-2021 23:10:13 (GMT)	Plant Review and Implementation
Workflow - LBL Approval with Proof Review Apr 7, 2021	07-Apr-2021 05:51:40 (GMT)	Milada Palkociova	Quality Approval	09-Apr-2021 09:25:01 (GMT)	Plant Review-1
Workflow - LBL Approval with Proof Review Apr 7, 2021	07-Apr-2021 05:51:40 (GMT)	Milada Palkociova	Quality Approval	09-Apr-2021 09:21:21 (GMT)	Proof Review

Workflow - LBL Approval with Proof Review Apr 7, 2021	07- Apr-2021 05:51:40 (GMT)	Milada Palkociova	Final Quality Approval	07- Apr-2021 08:29:42 (GMT)	Final Artwork Check
Workflow - LBL Approval with Proof Review Mar 31, 2021	31- Mar-2021 09:12:36 (GMT)	Milada Palkociova	Final Quality Approval	06- Apr-2021 13:12:00 (GMT)	Final Artwork Check