

# TECHNICAL DATA SHEET



## INCONTINENCE PADS

### SENI MAN



## I. PRODUCT DESCRIPTION

Seni Man incontinence pads for men with light and moderate incontinence. Anatomically shaped absorbent core with cellulose and superabsorbent beads. For better absorption and side leakages prevention Seni Man Normal and Extra pads are equipped with distributing layer [EDS] and side gathers with elastic yarn, and Seni Man Super pads have side folds with elastic yarn. Vapour permeable outer layer on the whole surface.

## II. RAW MATERIALS

Cellulose, superabsorbent, vapour permeable laminate, hydrophilic topsheet nonwoven, acquisition layer (not apply to the Super size), hydrophobic nonwoven (not apply to the Super size), elastic yarn, hot-melt adhesives, release paper, pouch film (not apply to the Super size).

## III. DIMENSIONS

| Size   | Total length* [mm] | Length of absorbent core* [mm] | Total width* [mm] |      | Szerokość wkładu chłonnego* [mm] |      |
|--------|--------------------|--------------------------------|-------------------|------|----------------------------------|------|
|        |                    |                                | front             | back | front                            | back |
| Normal | 285                | 255                            | 215               | 74   | 185                              | 47   |
| Extra  | 285                | 255                            | 215               | 74   | 185                              | 47   |
| Super  | 400                | 320                            | 90                | 90   | 90                               | 90   |

\* tolerance of dimensions: +/- 5 mm

## IV. PROPERTIES

| Size   | Weight [g] no less than | Absorbency [g] no less than (ISO 11948-1) |
|--------|-------------------------|-------------------------------------------|
| Normal | 15,2                    | 300                                       |
| Extra  | 19                      | 450                                       |
| Super  | 33,5                    | 820                                       |

## V. BASIC LOGISTIC DATA

| Size   | Code            | Package [pcs.] | Packages per box | Boxes per Pallet | Box type |
|--------|-----------------|----------------|------------------|------------------|----------|
| Normal | SE-095-MN15-G03 | 15             | 10               | 84               | carton   |
| Extra  | SE-095-EX15-G03 | 15             | 10               | 70               | carton   |
| Super  | SE-095-SU10-G03 | 10             | 12               | 42               | carton   |
|        | SE-095-SU20-G03 | 20             | 6                | 42               | carton   |

The incontinence pads must be placed on wooden pallet 1200x800 mm, covered by PE sheets (from the bottom and from the top) and wrapped together with PE stretch film.

## VI. RECOMMENDED STORAGE CONDITIONS


It is recommended to store the incontinence pads in a dry and well ventilated warehouse at stable temperature between 5°C and 35°C, at humidity between 25% and 75% and in their original packagings. The packagings must remain intact without openings in stretch film in order to avoid dust, odour and/or other contamination.

## VII. OTHER

Validity: 60 months.

**Issued by:** Piotr Przytarski, Kierownik Działu Rozwoju Technologii

**Approved by:** Jan Rosinski, Główny Technolog Grupy TZMO



Date of issue: 12.06.2020