

TECHNICAL DATA SHEET



HYGIENIC UNDERPADS

SENI SOFT SUPER



I. PRODUCT DESCRIPTION

Hygienic underpads are designed to be used as additional protection for people with incontinence. They allow for keeping bedding, bed, wheelchair, chair, etc. clean. Bottom outer layer is made of white, impermeable, anti-slip PE film (with printed pattern on). Absorbent core is made of cellulose fluff pulp and is additionally covered with hygienic tissue. Top layer is made of nonwoven.

II. RAW MATERIALS

Cellulose, hydrophilic nonwoven, PE film, hygienic tissue, hot- melt adhesives.

III. DIMENSIONS*

Size [cm]	Total length [mm]	Length of absorbent core [mm]	Total width [mm]	Width of absorbent core [mm]
40x60	400	320	600	580
60x60	600	520	600	580
90x60	900	820	600	580

* tolerance of dimensions +/- 10 mm

IV. PROPERTIES

Size [cm]	Weight [g] no less than	Absorbency [g] no less than (ISO 11948-1)
40x60	47,9	850
60x60	64,9	1400
90x60	95,0	2000

V. BASIC LOGISTIC DATA

Code	Size [cm]	Package [pcs]	Packages per box	Boxes per pallet	Box type
SE-091-SU05-001	40x60	5	24	24	carton
SE-091-SU30-001	40x60	30	4	36	carton
SE-091-SU05-002	60x60	5	24	20	carton
SE-091-SU05-SK2	60x60	5	24	20	carton

SE-091-SU30-002	60x60	30	4	24	carton
SE-091-SU30-SK2	60x60	30	4	24	carton
SE-091-SU05-003	90x60	5	12	28	carton
SE-091-SU05-SK3	90x60	5	12	28	carton
SE-091-SU05-I03	90x60	5	12	28	carton
SE-091-SU30-003	90x60	30	4	24	carton
SE-091-SU30-SK3	90x60	30	4	24	carton
SE-091-SU30-I03	90x60	30	4	24	carton
SE-091-SU25-D03	90x60	25	2	48	carton

The underpads must be placed on wooden pallet 1200x800 mm, covered by PE sheets (from the bottom and from the top) and wrapped together with PE stretch film.

VI. RECOMMENDED STORAGE CONDITIONS


It is recommended to store underpads in a dry and well ventilated warehouse at stable temperature between 5°C and 35°C, at humidity between 25% and 75% and in their original packagings. The packagings must remain intact without openings in stretch film in order to avoid dust, odour and/or other contamination.

VII. OTHER

Medical device marked with CE, class I.

Validity: 60 months.

Issued by: Marta Piątkowska, R&D Specialist



Approved by: Tomasz Przybylski, Deputy Director of Production and Innovation of TZMO Group

Date of issue: 19.04.2021