

4560-164 Guilhufe, Penafiel - Portugal

Ref: 485-557

## **Product Specification**

Brand:



## DRQ 14/21 REV.00 2020/09/08

Description:	Absorbent pad non-adherent

Final dimensions:	20 cm x20 cm	
Composition:	Cellulose, Cellulose pulp, Polyethylene, Polypropylene nonwoven, SSMMS blue PP NW (5 layers).	
Construction:	Abs.Pad non adherent	
Latex:	Free	
PVC:	Free	
Re use:	No	
Width variation:	± 5 mm	
Length variation:	+2/-0 mm	
Colour:	white/ blue	
Sterilization method:	Ethylene Oxide	
Shelf life:	60 months	
CFU/g	N.A.	



CFU/g	N.A.	D)
Nr and Revision Date:	17 2021/02/10	
CE Class:	I Steril	
CE mark:	CE1639	
Compliance with:	MDD 93/42/EEC amended by 2007/47/EC	
	ISO_10993-1.2018; EN_ISO_11607_1 2017 ASTM_F1980_16; ISO_14644-1:2015; ISO_	AC18;EN 556-1.2001_AC.2006;ISO_11135.2014; ;EN_ISO_11607-2.2017;ISO_TS_16775.2014; _14698_1_2003;EN_ISO_14971.2012;EN_ISO_15223-
Applicable norms:	1.2016; EN_1041.2008_A1.2013;	

Level		Туре	Dimension		Stamp	Labels
1	1	Flexible peel-open pack	134x240mm			
2	30	Folding box	142x315x250m	m		labeled
3	180	Carton box	780x285x315m	m		labeled

The device is intended to be used on a single patient during one single procedure to prevent risks from contamination, including possible infections, and/or functional damage.	
Store in a dry place at room temperature and protect from sunlight. The product has a good stability and	

Storage instructions	Store in a dry place at room temperature and protect from sunlight. The product has a good stability and therefore, normal conditions of medical devices storage shall apply.	
Waste instructions	No special requirements. The waste must be disposed according to local legal regulations for the disposable of medical waste.	
	Primary dressing intended for use on moderately to heavily exuding wounds and open drains. Universal use as absorbent pad in general patient care. Ideal for situations where a great capacity of absorption is needed, such as bleedings and chronicle wounds with moderate to high exudation. Secondary dressings for high exuding wounds.	
Intended use		
	Any modification of the product as presented or any use out of its intended use are not of the manufacturer responsibility.	

Non-harmonized symbols used in article labeling	non adherent	non-adherent
	MD	Medical Device

Environment - Packaging waste   Weight (gr) o	f packing material used in 100 pcs 1309
---	---

## ISSUED BY: Quality Management Department, 2021/05/14