


TECHNICAL DATA SHEET



BREATHABLE ALL-IN-ONE DIAPERS

SUPER SENI QUATRO



I. PRODUCT DESCRIPTION

All-in-one diapers for adults with very heavy urinary and faecal incontinence. Two anatomically shaped absorbent cores of cellulose pulp with SAP. Distributing layer (EDS). Side gathers along absorbent core are directed outside. Elastic leg cuffs preventing side leakages in groin area. Elastic waistband at the front and back. Four elastic, multi fastening combi-tapes (hook&loop fastener and adhesive tape). Wetness indicator: when diaper is exposed to wetness, ink inscription fades and two yellow adhesive stripes change colour into green. Vapour permeable outer layer on the whole surface.

II. RAW MATERIALS

Nontreated pulp fluff, superabsorbent, vapour permeable white laminate with claret stripes, topsheet hydrophilic-hydrophobic nonwoven, hydrophobic nonwoven, acquisition layer – white (EDS), hygienic tissue, elastic yarn, elastic closure tapes, elastic waistband, hot-melt wetness indicator, ink of indication, hot-melt, natural green tea extract.

III. DIMENSIONS

Size	Total length* [mm]	Length of absorbent core* [mm]		Total width* [mm]	Width of absorbent core* [mm] (front-crotch-back)	Suggested waist/hips girth [cm]
		top	bottom			
Small	720	370	580	520	250-150-250	55-80
Medium	830	400	690	640	260-160-260	75-110
Large	970	490	820	790	280-160-280	100-150
Extra Large	970	490	820	860	280-160-280	130-170

* tolerance of dimensions: +/- 10 mm

IV. PROPERTIES

Size	Weight [g] no less than	Absorbency [g] no less than (ISO 11948-1)	Retention [g] average (NAFC)
Small	99,9	2700	920
Medium	158,1	3700	1280
Large	170,7	4000	1400
Extra Large	172,6	4000	1400

V. BASIC LOGISTIC DATA

Size	Code	Package [pcs]	Packages per box	Boxes per pallet	Box type
Small	SE-094-SM10-A04	10	12	16	carton
Medium	SE-094-ME10-S04	10	6	24	carton
Large	SE-094-LA10-A04	10	6	24	carton
Extra Large	SE-094-XL10-A04	10	6	24	carton

The all-in-one diapers must be placed on wooden pallet 1200x800 mm, covered by PE sheet (from the bottom and from the top) wrapped together with PE stretch film.

VI. RECOMMENDED STORAGE CONDITIONS

It is recommended to store all-in-one diapers in a dry and well ventilated warehouse, at stable temperature between 5°C and 35°C and humidity level between 25% and 75%, in their original packaging. The packaging must remain intact without openings in stretch film in order to avoid dust and/or other contaminations.

VII. OTHER

Medical device marked with CE, class I.

Validity: 60 months

Issued by: Aleksandra Anklewicz, Technologist

A. Anklewicz

Approved by: Tomasz Przybylski, Deputy Director for
Production and Innovation of TZMO Group

Date of issue: 23.09.2020