

TECHNICAL DATA SHEET



DIAPERS WITH HIP BANDS

SENI OPTIMA PLUS



I. PRODUCT DESCRIPTION

Diapers with hip bands for adults with heavy urinary and faecal incontinence. Two anatomically shaped absorbent cores of cellulose pulp with SAP. Distributing layer (EDS). Side gathers along absorbent core are directed outside. Elastic leg cuffs preventing side leakages in groin area. Elastic waistband at the front and back. Two nonwoven belts, one of which has hook fastener. Two multi fastening combi-tapes (hook & loop fastener & adhesive tape). Wetness indicator: when diaper is exposed to wetness two yellow stripes change colour into green. Vapour permeable outer layer on the whole surface.

II. RAW MATERIALS

Nontreated pulp fluff, superabsorbent, vapour permeable white laminate with print, topsheet hydrophilic-hydrophobic nonwoven, hydrophobic nonwoven, acquisition layer – white (EDS), hygienic tissue, elastic yarn, closure tapes, hook fastener, belt nonwoven, elastic waistband, ink of indication, hot-melt wetness indicator, hot-melt, natural green tea extract.

III. DIMENSIONS

Size	Total length* [mm]	Length of absorbent core* [mm]		Total width* [mm]	Width of absorbent core* [mm] (front-crotch-back)	Suggested waist/hips girth [cm]
		top	bottom			
Medium	790	340	590	350	220-160-220	70-105
Large	840	400	690	415	245-160-245	80-125
Extra Large	990	420	740	520	300-170-340	105-155

* tolerance of dimensions: +/- 10 mm

IV. PROPERTIES

Size	Weight [g] no less than	Absorbency [g] no less than (ISO 11948-1)	Retention [g] average (NAFC)
Medium	118,7	2600	740
Large	159,6	2900	760
Extra Large	170,0	3500	1200

V. BASIC LOGISTIC DATA

Size	Code	Package [pcs]	Packages per box	Boxes per pallet	Box type
Medium	SE-097-ME10-G02	10	6	30	carton
Large	SE-097-LA10-G02	10	6	24	carton
Extra Large	SE-097-XL10-G02	10	8	16	carton

Diapers with hip bands must be placed on wooden pallet 1200x800 mm, covered by PE sheet (from the bottom and from the top) wrapped together with PE stretch film.

VI. RECOMMENDED STORAGE CONDITIONS

It's recommended to store diapers with hip bands for adults in a dry and well ventilated warehouse, at stable temperature between 5°C and 35°C and humidity level between 25% and 75%, in their original packaging. The packaging must remain intact without openings in stretch film in order to avoid dust and/or other contaminations.

VII. OTHER

Medical device marked with CE, class I.

Validity: 60 months

Issued by: Agnieszka Tarka, R&D Specialist



Approved by: Tomasz Przybylski, Deputy Director for
Production and Innovation of TZMO Group

Date of issue: 29.09.2020