


Portex™ Single-Use Tracheal Tube Exchange Guides

Sizes 5CH and 10CH



Straight tip

Device Classification

GMDN: 42084

Classification: Is

Device Description

The Portex single-use tracheal tube exchange guides are single-use, sterile tracheal tube guides designed for securing access to the patient's airway during changing of a tracheal tube. The tracheal tube guides are intended for use with the following tracheal tube sizes:

Tracheal Tube Guide (Ch)	5	10
Tracheal Tube Size (mm) I.D.	2.5–3.5	4.0–6.0

The depth of insertion of the guide can be easily determined by using the graduation marks at 10, 20 and 30 cm from the guide tip.

Precautions

1. The guide should not be used for nasal intubation.
2. The tracheal tube should not be threaded over the guide into the trachea without a laryngoscope in place.

Component Composition

- > 5 Ch tube: Polythene HD
- > 10 Ch tube: Polythene HD/ Yellow Polyurethane

This device is not made with plasticiser diethylhexylphthalate (DEHP).

This device is not made with natural rubber latex.

This device is not manufactured using derivatives of tissues or cells of animal origin.

Legal Manufacturer Name and Address

ICU Medical, Inc.

6000 Nathan Lane North

Minneapolis, MN 55442 USA

Country of Origin

Mexico

Sterilisation Method

Ethylene Oxide (EO)

Shelf Life

5 years

Labelling and Packaging

Container Type	Length (mm)	Width (mm)	Height (mm)
100/125/005			
Single unit pouch	725	70	
Shelf Carton (10 units)	725	68	40
Transit Carton (200 units)	731	352	168
100/125/010			
Single unit pouch	840	70	
Shelf Carton (10 units)	836	70	42
Transit Carton (100 units)	845	225	146

One unit packaged per pouch. Ten units packaged per shelf carton. For 100/125/005, twenty shelf cartons packaged per transit carton, for a total of two-hundred units. For 100/125/010, ten shelf cartons packaged per transit carton, for a total of one-hundred units. Discard if open, wet or damaged. Product should be kept out of direct sunlight.

The lot number, manufacturing date and expiration date are located on the pouch, shelf carton and transit carton labels.

Product Specifications

Portex Single-Use Tracheal Tube Exchange Guides

Product Code	Recommended Tube I.D. (mm)	Size Ch	O.D. (mm)	Length (mm)	Units per case
100/125/005	2.5–3.5	5	1.7	500	10
100/125/010	4.0–6.0	10	3.3	700	10

The product complies with current legislation and has the corresponding CE marking (2797).
For additional information, warnings and/or safety precautions, refer to the manufacturer's Instructions for Use.