

TECHNICAL DATA SHEET



SENI LADY SLIM bladder control liners

SENI LADY bladder control pads



I. PRODUCT DESCRIPTION

Seni Lady Slim bladder control liners/Seni Lady bladder control pads designed for adults with very light, light and moderate incontinence. Anatomically shaped absorbent core with cellulose and superabsorbent beads. Distribution layer (EDS) (except Micro size) and embossing support urine distribution. Purple side gathers with elastic yarn along absorbent core are directed outside prevent side leakages (except Micro size). Vapour permeable outer layer covering entire surface.

II. RAW MATERIALS

Cellulose, superabsorbent, vapour-permeable laminate, hydrophilic nonwoven, hydrophobic nonwoven (except Micro size), distribution layer (EDS) (except Micro size), flat elastic (except Micro, Mini, Mini Plus, Normal sizes), elastic yarn (except Micro size), hot-melt adhesives, release paper, water-soluble paint.

III. DIMENSIONS

Size	Total length [mm]*	Length of absorbent core [mm]*	Total width [mm]** (front – crotch – back)	Width of absorbent core [mm]** (front – crotch – back)
Slim Micro	185	167	72-59-72	58-45-58
Slim Mini	225	205	95-80-95	75-60-75
Slim Mini Plus	255	225	100-85-100	80-65-80
Slim Normal	280	245	105-90-105	85-70-85
Slim Extra	295	265	130-110-140	100-75-110
Extra Plus	330	295	140-115-150	110-80-115
Super	370	335	165-130-180	130-90-145
Plus	420	380	185-145-205	145-105-165

* tolerance of dimensions: +/- 3 mm (Micro); +/- 5 mm (Mini, Mini Plus, Normal, Extra, Extra Plus, Super, Plus).

** tolerance of dimensions: +/- 1 mm (Micro); +/- 5 mm (Mini, Mini Plus, Normal, Extra, Extra Plus, Super, Plus).

IV. PROPERTIES

Size	Weight [g] no less than	Absorbency [g] no less than (ISO 11948-1)	Retention [g] average (NAFC)
Slim Micro	2,95	60	30
Slim Mini	8,3	200	70
Slim Mini Plus	10,1	260	115
Slim Normal	12,3	340	150
Slim Extra	15,9	470	200
Extra Plus	22,5	590	220
Super	32,2	800	320
Plus	37,0	950	370

V. BASIC LOGISTIC DATA

Size	Code	Package [pcs]	Packages per box	Boxes per pallet	Box type
Slim Micro	SE-095-MC20-DE1	20	16	128	carton
Slim Mini	SE-095-MI20-DE1	20	30	36	carton
Slim Mini Plus	SE-095-MP20-DE1	20	30	30	carton
Slim Normal	SE-095-NO20-DE1	20	18	45	carton
Slim Extra	SE-095-EX15-DE1	15	24	30	carton
Extra Plus	SE-095-EP15-DE1	15	20	24	carton
Super	SE-095-SU15-DE1	15	16	24	carton
Plus	SE-095-PL15-DE1	15	16	24	carton

The bladder control liners/bladder control pads must be placed on wooden pallet, covered by PE sheets (from the bottom and from the top) and wrapped together with PE stretch film.

VI. RECOMMENDED STORAGE CONDITIONS

It is recommended to store bladder control liners/bladder control pads in a dry and well ventilated warehouse at stable temperature between 5°C and 35°C, at humidity between 25% and 75% and in their original packagings. The packagings must remain intact without openings in stretch film in order to avoid dust, odour and/or other contamination.

VII. OTHER

Medical device marked with CE, class I.

Validity: 60 months.

Issued by: Marta Stachlewska, R&D Specialist

Marta Stachlewska
Marta Stachlewska
Technolog

Approved by: Tomasz Przybylski, Deputy Director of
Production and Innovation of TZMO Group

Date of issue: 19.04.2021