

# SOL-M™ Catheter Tip Syringe without Needle



## SOL-M™ Catheter Tip Syringe without Needle is used to:

- irrigate
- withdraw fluid from a body cavity or wound
- instil fluid into a body cavity or wound


## Features:

- **High-transparency barrel** - makes it easy to see syringe content.
- **Silicone lubrication** - ensures smooth plunger movement.
- **Safe plunger backstop** - prevents accidental withdrawal of plunger.
- **Large and stable finger grip** - ensures stability and comfort during injection.
- **Bold scale calibrated in oz and ml** - to ensure precise dosage. Allows multiple technique and practices.




**SOL-MILLENNIUM®**  
Building a Healthier Tomorrow

# SOL-M™ Catheter Tip Syringe without Needle



## Tip cover cap for 60 ml

Protect the syringe content during transportation from filling site to fluid discharge site.



## Luer adapter for 100ml

Allows filling the syringe from sterile solution using a standard Luer connection.

- Not made with natural rubber latex.
- DEHP-, PVC- free.
- Sterile and single-use.
- Sterilization method: ethylene oxid.
- Blister packaging.
- Unique Device Identification (UDI) on the blister.

Reference number	Size	Shelf Box/Shipping Case (pcs)
180060CT	60 ml	30/240
180100CT	100 ml	30/60

CE<sub>2797</sub>

### For more information, please contact:

Sol-Millennium Europe Sp. z o.o  
Twarda 18, 00-105 Warsaw, Poland  
Tel: +48 22 202 67 23  
[info.EU@sol-m.com](mailto:info.EU@sol-m.com)  
[sol-m.com](http://sol-m.com)

EMEA-MD-231-0821

 **SOL-MILLENNIUM®**  
Building a Healthier Tomorrow