

TECHNICAL DATA SHEET



ELASTIC DISPOSABLE UNDERWEAR

SENI ACTIVE TRIO



I. PRODUCT DESCRIPTION

Seni Active Trio elastic disposable underwear are designed for people with heavy urinary incontinence. Top and bottom layers are made of soft non-woven. Absorbent core is made of cellulose pulp combined with superabsorbent beads. Contains distribution layer (EDS). Anti-leakage protection is provided by side standing gathers along absorbent core and leg cuffs with elastic yarn in groin area. Wetness indicator – black, printed inscription that fades when exposed to wetness. Vapour permeable outer layer on the whole surface. Z-fold tape in back area to be wrapped around the underwear after use.

II. RAW MATERIALS

Cellulose, hydrophilic nonwoven, hydrophobic nonwoven, distribution layer (EDS), vapour permeable film, elastic yarn, hot-melt adhesives, superabsorbent, Z-fold tape.

III. DIMENSIONS*

Size	Total length [mm]	Length of absorbent core [mm]	Total width [mm]	Width of absorbent core [mm] (front – crotch – back)	Suggested waist/hips girth [cm]
Medium	810	510	660	160-100-160	80-110
Large	870	510	770	160-100-160	100-135
Extra Large	920	510	890	160-100-160	120-160

* tolerance of dimensions: +/- 10 mm

IV. PROPERTIES

Size	Weight [g] no less than	Absorbency [g] no less than (ISO 11948-1)	Retention [g] average (NAFC)
Medium	95,5	2500	1100
Large	97,8	2500	1100
Extra Large	103,5	2500	1100

V. BASIC LOGISTIC DATA

Size	Code	Package [pcs]	Packages per box	Boxes per pallet	Box type
Medium	SE-096-ME10-A03	10	8	24	carton
Large	SE-096-LA10-A03	10	8	24	carton
Extra Large	SE-096-XL10-A03	10	6	28	carton

The disposable underwear must be placed on wooden pallet, covered by PE sheets from the top and PE sheets or cardboard divider sheet from the bottom and wrapped together with PE stretch film.

VI. RECOMMENDED STORAGE CONDITIONS

It is recommended to store disposable underwear in a dry and well ventilated warehouse at stable temperature between 5°C and 35°C, at humidity $\leq 70\%$ and in their original packagings. The packagings must remain intact without openings in stretch film in order to avoid dust, odour and/or other contamination.

VII. OTHER

Medical device marked with CE, class I.

Validity: 60 months.

Issued by: Justyna Dziwulska, Technology Specialist

Approved by: Tomasz Przybylski, Authorized Proxy/
Deputy Director of Production and Innovation

Date of issue: 11.08.2022

this is legal evidence of all created e-signatures and e-seals

Certificate for a document with an Autenti ID: 22128df4-ba51-4f93-a56f-74d8e92100c9
created at: 2022-08-11 08:16 (GMT+02:00)