

TECHNICAL DATA SHEET



DIAPERS WITH HIP BANDS

SENI OPTIMA QUATRO



I. PRODUCT DESCRIPTION

Diapers with hip bands for adults with severe urinary and faecal incontinence. Two anatomically shaped absorbent cores of cellulose pulp with SAP. Distributing layer (EDS). Side gathers along absorbent core are directed outside. Elastic leg cuffs preventing side leakages in groin area. Elastic waistband at the front and back. Two nonwoven belts, one of which has hook fastener. Two multi fastening combi-tapes (hook & loop fastener & adhesive tape). Wetness indicator: when diaper is exposed to wetness two yellow stripes change colour into green. Vapour permeable outer layer on the whole surface.

II. RAW MATERIALS

Nontreated pulp fluff, superabsorbent, vapour permeable white laminate with print, topsheet hydrophilic-hydrophobic nonwoven, hydrophobic nonwoven, acquisition layer – white (EDS), elastic yarn, closure tapes, hook fastener, belt nonwoven, elastic waistband, ink of indication, hot-melt wetness indicator, hot-melt, natural green tea extract.

III. DIMENSIONS

Size	Total length* [mm]	Length of absorbent core* [mm]		Total width* [mm]	Width of absorbent core* [mm] (front-crotch-back)	Suggested waist/hips girth [cm]
		top	bottom			
Medium	790	340	590	350	220-160-220	70-105
Large	840	400	690	415	245-160-245	80-125

* tolerance of dimensions: +/- 10 mm

IV. PROPERTIES

Size	Weight [g] no less than	Absorbency [g] no less than (ISO 11948-1)	Retention [g] average (NAFC)
Medium	140,6	3500	1470
Large	189,0	4300	1600

V. BASIC LOGISTIC DATA

Size	Code	Package [pcs]	Packages per box	Boxes per pallet	Box type
Medium	SE-097-ME10-G04	10	6	30	carton
Large	SE-097-LA10-G04	10	6	24	carton

Diapers with hip bands must be placed on wooden pallet 1200x800 mm, covered by PE sheet (from the bottom and from the top) wrapped together with PE stretch film.

VI. RECOMMENDED STORAGE CONDITIONS

It is recommended to store diapers with hip bands for adults in a dry and well ventilated warehouse, at stable temperature 5°C - 35°C and humidity level ≤ 70%, in their original packaging. The packaging must remain intact without openings in stretch film in order to avoid dust, foreign odours and/or other contamination.

VII. OTHER

Medical device marked with CE, class I.

Validity: 60 months

Issued by: Aleksandra Anklewicz
Technologist

Approved by: Tomasz Przybylski
Authorized Proxy / Deputy Director of
Production and Innovation

Date of issue: 15.06.2022

this is legal evidence of all created e-signatures and e-seals

Certificate for a document with an Autenti ID: 105bc8bf-d3d2-4490-b8dc-2324a6045fd7
created at: 2022-06-15 14:24 (GMT+02:00)

