


PRODUCT TECHNICAL DATASHEET

ITEM REF:	011-42394-12	REVISION:	01
-----------	--------------	-----------	----

LEGAL MANUFACTURER	ICU Medical, Inc. 951 Calle Amanecer, San Clemente, CA 92673, USA						
IMPORTER	ICU Medical B.V. Hofspoor 3, 3994 VZ Houten, The Netherlands						
EU AUTHORISED REPRESENTATIVE (EUAR)	Medical Device Safety Service GmbH (MDSS), Schiffgraben 41, 30175 Hannover, Germany						
ASSEMBLY SITE	ICU Medical de Mexico, S.A. de C., Avenida Cuarza No. 250, Colonia Rancho Santa Clara, Delegacion Manadero Ensenada, Baja California, Mexico 22790						
CLASSIFICATION CODE	GMDN Code: 35504						
INTENDED USE	Blood pressure transducer holder. A device used to secure a pressure transducer used for invasive blood pressure measurement at a desired height. The holder can be fastened to (e.g., a bed or an operating table) and can be adjusted so that the transducer is at the same level as the heart of the patient. This is a reusable device.						
ITEM DESCRIPTION	MOUNT HOLDER FOR POLE MOUNTING MONITORING KIT COMPONENTS Reusable Device: Do Not Discard						
PRIMING VOLUME (ml)	n/a	LENGTH (cm)	23.49	WEIGHT (g)	190.40	CASE QTY	5
							
LIST OF COMPONENTS	1	MANIFOLD HOLDER MH02	STAINLESS STEEL (NON-FLUID PATH) PVC (NON-FLUID PATH)				
STERILIZATION AND SHELF LIFE	NON-STERILE						
MATERIAL COMPLIANCE AND BIOCOMPATIBILITY	This product is made of non-latex, and non-DEHP components.						
TRACEABILITY	Lot number provides full traceability of all components and manufacturing processes.						
STORAGE	Store in a dry and clean place. Product should be retained in packaging until ready for use.						
DISPOSAL	The user must dispose of the device according to hospital disposal policy.						