



Zelletten® Cellulose swabs off the roll

REF 13349, 13356

January 2023

1. Composition of the product

Zelletten® Cellulose Swabs of the roll consist of:

- bleached, 12-sheet cellulose wadding

This product data sheet covers the following items:

REF	Product description	Feature	Package size
13349	Zelletten® Cellulose Swabs	5 x 4 cm	Roll with 300 pieces in transparent dispenser
13356		5 x 4 cm	2 rolls á 500 Cellulose Swabs

2. Packaging, structure and composition

2.1 Unit container

see 2.2

2.2 Shelf container

REF 13349

1 Roll in a transparent dispenser made of polysterol
folding box, made of cellulose
Packaging units see product table above

REF 13356

bag made of polyethylene with drawcord-clip
Packaging units see product table above

2.3 Transit container

1 corrugated cardboard box (cellulose)

REF	packings	units
13349	20	20 folding boxes
13356	12	12 bags = 24 rolls

3. Manufacture

Zelletten® Cellulose Swabs of the roll are produced according to specification under hygienic conditions and wrapped as described in the relevant packaging specification.

4. Description

White, soft cellulose swabs with typical cellulose odour.

5. Properties

Zelletten® Cellulose Swabs of the roll are bleached without chlorine.

6. Intended purpose

Zelletten® Cellulose Swabs are used to apply ointments and to swab body fluids when treating minor injuries.

7. Medical device classification

Zelletten® Cellulose Swabs are medical devices of class I according to rule 4 Indent 1. (acc. to Regulation (EU) 2017/745)

8. Biological evaluation and biocompatibility (DIN EN ISO 10993)

The starting materials used in the manufacture of Zelletten® Cellulose Swabs of the roll are safe if the product is used appropriately.

The purpose of this documentation and the statements made therein is to show that there is no risk involved in the use of the medical device Zelletten® Cellulose Swabs and that it is designed, manufactured and packaged in such a way that it will not compromise the clinical condition or the safety of patients, or the safety and health of users and other persons when used under the conditions and for the purposes intended.

9. Requirements regarding the safety/harmlessness of the constituents

According to the supplier's information, there are no ingredients of concern.

10. Stability

Stored appropriately (keep dry, protected from dust and sun) Zelletten® Cellulose Swabs off the roll have a shelf life of 5 years from the date of manufacture.

11. Disposal

The user is advised to observe current national legislation, norms and guidelines, regulating the disposal of medical refuse.
Packaging materials must also be disposed of in compliance with applicable national requirements.

Lohmann & Rauscher International GmbH & Co. KG
D-56579 Rengsdorf
signed by

Peter Grillitsch
Director Global Marketing