

TECHNICAL DATA SHEET



ALL-IN-ONE DIAPERS FOR ADULTS

SENI SUPER PLUS XXL



I. PRODUCT DESCRIPTION

All-in-one diapers for adults with heavy urinary and faecal incontinence.

Two anatomically shaped absorbent cores of cellulose pulp with SAP.

Distributing layer (EDS).

Double protection against leaks : inner side gathers along absorbent core are directed outside; outside elastic leg cuffs in groin area.

Elastic waistband at the front and back.

Four multi fastening adhesive tapes.

Two nonwoven belts on the extension of the side wings.

Vapour permeable outer layer on the whole surface.

Wetness indicator: when diaper is exposed to wetness two yellow adhesive stripes change color into green-blue.

The purple print on the insulating layer has an informative function :

- product name and size,
- 7 drops and purple color – determination of the absorbency level dedicated to heavy urinary,
- pictograms that inform how to put the diaper: A – standing, B – lying down,
- transverse lines dividing the indicator into zones making it easier to assess the degree of diaper filling.

Additional ink of print defining date and time of product production.

II. RAW MATERIALS

Nontreated pulp fluff, superabsorbent, vapour permeable white laminate with purple print, topsheet hydrophilic-hydrophobic nonwoven, hydrophobic nonwoven, acquisition layer – white (EDS), hygienic tissue, elastic yarn, tapes, belt nonwoven, elastic waistband, hot-melt wetness indicator, ink of indication, hot-melts, natural green tea extract.

III. DIMENSIONS

Size	Total length* [mm]	Length of absorbent core* [mm]		Total width* [mm] (front*-back w/ belt)**	Width of absorbent core* [mm] (front-crotch-back)	Suggested waist/hips girth [cm]
		top	bottom			

XXL	1015	530	870	640-1620	280-160-280	160-210
-----	------	-----	-----	----------	-------------	---------

* tolerance of dimensions: +/- 10 mm

** tolerance of dimensions: +/- 20 mm

IV. PROPERTIES

Size	Weight [g] no less than	Absorbency [g] no less than (ISO 11948-1)	Retention [g] average (NAFC)
XXL	152	3200	1150

V. BASIC LOGISTIC DATA

Size	Code	Package [pcs]	Packages per box	Boxes per pallet	Box type
XXL	SE-094-2X10-G02	10	6	24	carton

The all-in-one diapers must be placed on wooden pallet 1200x800 mm, covered by PE sheet (from the bottom and from the top) wrapped together with PE stretch film.

VI. RECOMMENDED STORAGE CONDITIONS

It is recommended to store all-in-one diapers in a dry and well ventilated warehouse, at stable temperature 5°C - 35°C and humidity level ≤ 70%, in their original packaging. The packaging must remain intact without openings in stretch film in order to avoid dust, foreign odours and/or other contamination.

VII. OTHER

Medical device marked with CE, class I.

Validity: 60 months

Issued by: Agnieszka Tarka
R&D Specialist

Approved by : Magdalena Ćwiertnia
Leader of Product Technology for Incontinence
TZMO SA

Date of issue: 12.07.2024

poświadczenie złożenia podpisów i pieczęci elektronicznych

Certyfikat dla dokumentu o Autenti ID: 3c367e75-4893-4e39-b678-a03af2776a5d
utworzonego: 2024-07-16 11:43 (GMT+02:00)